

TOPIC 3:

Sequence data


Outline

- 1. Different methods to acquire sequence data**
2. Understand sequence file formats
3. Preparing files for analysis
 - Tutorial looking at sequence data files and quality

Whole Genome Sequencing

Randomly sheer DNA and sequence all fragments

May use double-stranded nuclease treatment to reduce repetitive elements




Screen shot from the Integrated Genomics Viewer

Whole Genome Sequencing

Pros	Cons
All sites possible	Comparatively expensive per sample
Simple library prep	Storage and bioinformatics challenging with lots of samples

Pool Seq



Adapted from Fuentes-Pardo & Ruzzante 2017 Mol. Ecol

Pool Seq

Pros

All sites possible

Simple library prep

Cheaper than
individual WGS


Cons

Limited analysis
options

No haplotype
information

Best in cases where
 $\# \text{ samples} > \# \text{ reads}$

RNAseq



RNAseq

Pros

Many sites and only
in genes

Also get expression
information

Relatively easy to
assemble

Cons


Expression
differences
complicate SNP
calling

Expensive for pop
gen level sampling

Difficult library prep
(or so I'm told!)

Amplicon Sequencing

- Use PCR to amplify target DNA. Sequence many barcoded samples in one lane.
- Used to characterise microbiome by sequencing 16s rRNA



Amplicon Sequencing

Pros

Get incredible depth at single locus

Simple bioinformatics.

Cons

Limited to one or few loci


Mutations in primer site don't sequence

GT-seq

- Genotyping by Thousands
- Based on Amplicon sequencing
- Multiplex PCR amplify ~200 known SNPs and then sequence pooled PCR products.
- Very cheap (\$1/sample), and bioinformatically simple.
- Useful for genotyping thousands or tens of thousands of samples.
- Complicated initial set-up.

Sequence Capture

- Design probe sequences from genome resources, synthesis attached to beads
- Make WGS library, hybridize with probe set. Matching sequence will be captured, all others washed away
- Collect capture sequence, amplify and sequence



Sequence Capture

Pros

Relatively cheap
per sample

Good depth at
targeted sites

Cons

Requires
designing probes


Long library prep

Reduced Representation Sequencing

Instead of sequencing the whole genome, it can be sufficient to sequence just a part of it

A

RAD sequencing



B

double digest RADseq




Figure from Peterson et al PLoS One 2012

Reduced Representation Sequencing

Pros

Quick library prep for hundreds of samples


Comparatively cheap per sample cost

Cons

Relatively sparse SNPs compared to other methods
- limiting analysis options

Can have problems overlapping different library preps

Synthetic long reads



Barcodes read originating from individual DNA molecules


Sequence with Illumina reads

Original molecule can be reconstructed using the barcodes

Potentially very useful for genome assembly and phasing

Synthetic long reads - i.e. Linked Reads


E.g. Haplotagging



Bisulphite Sequencing

Unmethylated cytosines are converted to **U**racil

Methylated **C**pG sites are unchanged and are detected as polymorphisms



How to choose?

For example:

If you wanted to estimate demographic history from the distribution of allele frequencies, a reduced representation method might suffice to obtain an estimate of the site frequency spectrum


Or, if you want to perform a genome scan, looking at how haplotype frequencies varied among populations, you'd probably need deeper, whole genome information - it all depends on the questions you are tackling

Outline

1. Different methods to acquire sequence data
- 2. Understand sequence file formats**
3. Preparing files for analysis
 - Tutorial looking at sequence data files and quality

Part 2: Sequence file formats

Illumina sequencing




*Reverse strands are cleaved
after cluster formation*

*4 cycles are shown, but modern
Illumina machines are capable of
600 cycles in one run*

Part 2: Sequence file formats

Illumina sequencing




This process has generated 4 Images

Reverse strands are cleaved after cluster formation

4 cycles are shown, but modern Illumina machines are capable of 600 cycles in one run

Part 2: Sequence file formats

Illumina sequencing




In this example, there are:

8 Lanes, each containing
2 columns, within which there are
50 tiles, each tile is imaged 4 times/cycle
(4-channel chemistry)

Part 2: Sequence file formats

Illumina sequencing




So there are approximately $8 \times 2 \times 50 \times 4 = 3,000$ images generated per cycle

Each image is about 3Mb in size

For an Illumina run using 300 cycles, that would be $3000 \times 3 \times 300 = 2,700,000$ Mb of data (~2.7 Tb)

Part 2: Sequence file formats

Illumina sequencing



The number of channels refers to the numbers of colours the images detect

4-channel was Illumina's standard chemistry, but now 2-channel is more common

Part 2: Sequence file formats

Illumina sequencing

Using the stack of images from an Illumina machine you do the following:



1. Evaluates the light signal from every cluster to calculate the Quality Predictor Value (QPV), measuring things like:

The signal-to-noise ratio
Light Intensity


2. QPVs are converted into Phred quality scores (Q-scores) using a calibration curve built using previously sequenced samples

3. Convert the base call and the quality score into a FASTQ file

Part 2: Sequence file formats

Illumina sequencing

Using the stack of images from an Illumina machine you do the following:



You'll probably never do these steps yourself, but it's good to know where the data come from!

1. Evaluate the light signal from every cluster to

measuring

signal-to-noise ratio

Light Intensity

2. QPVs are converted into Phred quality scores (Q-scores) using a calibration curve built using previously sequenced samples

3. Convert the base call and the quality score into a FASTQ file

Part 2: Sequence file formats

Remember this from yesterday?

What a FASTA file looks like:

Sequence
name

```
>chr_1
TGGGCAAGGCTGATGAACAGCAGCTGCATAAATTCTCCCCTAATTATATTGTAATAGCT
GCAGCACAACAATAAAGCTTTGTTAGAGACATCTAGAGAATCACACACTGCATCTGTTCT
GCCGCTCTCCCTCTTGCTCTGTTCTGAGAAGCACTTGTTCACTGATTCTGGGTTTGTATT
TGTGTTTTTCATGCTTAACATTGTTATTTGTTTGCCTAGAAAGTTCTTTGATTGGGCCAA
ATTAGTCGATTTTAAAGAGTGCACTTCTCTAGTGCATGTAATCTATGTGGACATCTCAAT
AGCTGCTTAATTTGTTTAGTGGTAATCTCCTCTGAACAGAGAGAAAGGCCTACATGCAGC
CCTCAGAGGAGAGGTGTCAATCTCTCTTTGATTATCTCTTTGTTTCCCTTCAGAAGAATC
ATTCTAATCTGGTATTGTACAAGAGGAAATAAATGGGACTAAAACCAGGCATGCACCATC
TGATAGATTCACATCCCTAGAAGACTTTTGTTGTGTTTGTTCAGTGAGAGCCTGCTG
```

Nucleotide
sequence

FASTAs are plain text files

Part 2: Sequence file formats

Anatomy of a FASTQ file:

4 Lines instead of 2

```
@SRR357068.1 D042KACXX:3:1101:2690:2160 length=101
NCATCGTCCGGTATGTAGAACAGGGGAACCGGACGTTTTCCAAGGCGTAGC
CATGTTAGACAAGGCGCAGATATAGGTGA
+SRR357068.1 D042KACXX:3:1101:2690:2160 length=101
#4=DBDDDFHFFHHIGHIIJJJJJJJJJJJBHDAGHJGGGHIJHFFFFD
DEDCCDCCCCDDDDDBDBD>CDEE>C@CD
```

FASTQs are plain text files

Part 2: Sequence file formats

Anatomy of a FASTQ file:

1. Sequence
ID
(begins with “@” not “>”)

Typically contains information on the origin of the read - like which lane and tile it came from, where in the tile the cluster was located



```
@SRR357068.1 D042KACXX:3:1101:2690:2160 length=101
NCATCGTCCGGTATGTAGAACAGGGGAACCGGACGTTTTCCAAGGCGTAGC
CATGTTAGACAAGGCGCAGATATAGGTGA
+SRR357068.1 D042KACXX:3:1101:2690:2160 length=101
#4=DBDDDFHFFHHIGHIIJJJJJJJJJJJBHDAGHJGGGHIJHFFFFD
DEDCCDCCCCDDDDDBDBD>CDEE>C@CD
```

FASTQs are plain text files

Part 2: Sequence file formats

Anatomy of a FASTQ file:

1. Sequence
ID

(begins with “@” not “>”)

2. Nucleotide
sequence



```
@SRR357068.1 D042KACXX:3:1101:2690:2160 length=101
NCATCGTCCGGTATGTAGAACAGGGGAACCGGACGTTTTCCAAGGCGTAGC
CATGTTAGACAAGGCGCAGATATAGGTGA
+SRR357068.1 D042KACXX:3:1101:2690:2160 length=101
#4=DBDDDFHFFHHIGHIIJJJJJJJJJJJJJBHDAGHJGGGHIJHFFFFD
DEDCCDCCCCDDDDDBDBD>CDEE>C@CD
```

FASTQs are plain text files

Part 2: Sequence file formats

Anatomy of a FASTQ file:



FASTQs are plain text files

Part 2: Sequence file formats


Anatomy of a FASTQ file:



FASTQs are plain text files

Part 2: Sequence file formats

Anatomy of a FASTQ file:



FASTQs are plain text files

Part 2: Sequence file formats

Base quality scores (Q-scores)

$$Q_{Sanger} = -10\log_{10}(p)$$

Where p is the probability that a base call is incorrect

$$Q_{Solexa} = -10\log_{10}\left(\frac{p}{1-p}\right)$$

Remember, those probabilities are calculated using the QPVs in Illumina sequencing

Part 2: Sequence file formats


Base quality scores (Q-scores)

$$Q_{Sanger} = -10\log_{10}(p)$$

Where p is the probability that a base call is incorrect

$$Q_{Solexa} = -10\log_{10}\left(\frac{p}{1-p}\right)$$

Remember, those probabilities are calculated using the QPVs in Illumina sequencing



Red line is Sanger
Black line is Solexa

Asymptotically identical when p is small

Part 2: Sequence file formats

Base quality scores (Q-scores)

$$Q_{Sanger} = -10 \log_{10}(p)$$

What's the probability that the base is incorrect if Q=30?

Part 2: Sequence file formats

Base quality scores (Q-scores)

$$Q_{Sanger} = -10 \log_{10}(p)$$

What's the probability that the base is incorrect if Q=30?

$$p[Q30] = 0.001$$

$$p[Q20] = 0.01$$

$$p[Q10] = 0.1$$

Part 2: Sequence file formats

Base quality scores (Q-scores)

You probably noticed that the Q-scores in the FastQ files are not numeric


```
@SRR357068.1 D042KACXX:3:1101:2690:2160 length=101
NCATCGTCCGGTATGTAGAACAGGGGAACCGGACGTTTTCCAAGGCGTAGC
CATGTTAGACAAGGCGCAGATATAGGTGA
+SRR357068.1 D042KACXX:3:1101:2690:2160 length=101
#4=DBDDDFHFFHHIGHIIJJJJJJJJJJJJJBHDAGHJGGGHIJHFFFFD
DEDCCDCCCCDDDDDBDBD>CDEE>C@CD
```

Under Illumina sequencing,
ASCII encoding is used to refer to Q scores from 0 to 62

Slightly different encoding strategies are used by the different technologies

Part 2: Sequence file formats

PacBio Sequencing



Nanopore Sequencing




Figure 3. Single Molecule Sequencing Platforms

(A) Pacific Bioscience's SMRT sequencing. A single polymerase is positioned at the bottom of a ZMW. Phosphate-labeled versions of all four nucleotides are present, allowing continuous polymerization of a DNA template. Base incorporation increases the residence time of the nucleotide in the ZMW, resulting in a detectable fluorescent signal that is captured in a video.

(B) Oxford Nanopore's sequencing strategy. DNA templates are ligated with two adaptors. The first adaptor is bound with a motor enzyme as well as a tether whereas the second adaptor is a hairpin oligo that is bound by the HP motor protein. Changes in current that are induced as the nucleotides pass through the pore are used to discriminate bases. The library design allows sequencing of both strands of DNA from a single molecule (two-direction reads).

Outline

1. Different methods to acquire sequence data
2. Understand sequence file formats
- 3. Preparing files for analysis**
 - Tutorial looking at sequence data files and quality

Part 3: Preparing files for analysis

What do you do when you get your data?


Part 3: Preparing files for analysis

1) Check files for completeness, use checksums if file corruption is suspected

Downloading large data files takes a long time

There is a possibility of data corruption when files are downloaded

There are command line tools for verifying data complete-ness



MD5 and SHA-1 are the most common checksum methods

There is a short demonstration using SHA-1 sums in the tutorial

Part 3: Preparing files for analysis

- 1) Check files for completeness, use checksums if file corruption is suspected
- 2) Inspect quality statistics


There are many possible statistics to query:

- Number and length of sequences
- Base qualities**
- Poly A/T tails
- Presence of tag sequences (things that you added during library prep.)
- Sequence complexity (e.g. identify repetitive data ATATATATATATATATATA)


There are standard tools for examining these, such as prinseq and fastqc

Part 3: Preparing files for analysis

Distribution of base frequencies in GBS reads - with enzyme cut site




Distribution in RNAseq data - with no adapters/tags used



Part 3: Preparing files for analysis

A typical quality score distribution for Illumina reads



Part 3: Preparing files for analysis

1) Check files for completeness, use checksums if file corruption is suspected

2) Inspect quality statistics

3) Possible steps to clean files

- De-multiplex
- Trim adapters
- Filter/trim low quality base calls
- Remove duplicate sequences
- Remove contaminant sequences
- Remove sequences that are mainly adapter

Often done
by the
sequencing
centre

Important for
genotyping
and RNAseq

Important for
reference
assembly

Many programs to implement these steps!

Part 3: Preparing files for analysis

Quality trimming

Choice of quality score to filter to depends upon the application:

- Too low a quality score cutoff:
 1. increase run times and RAM usage
 2. Bad results (e.g. false SNP calls)
- Too high a quality score cutoff:
 1. Faster run times
 2. Potentially lose useful data (e.g. more fragmented assemblies or missing SNPs)


Q20 is a rule of thumb, but it depends on what you're doing

Part 3: Preparing files for analysis

Quality trimming

Del Fabbro et al 2013

Number of variants detected



Part 3: Preparing files for analysis

Contamination


Contamination in your samples can lead to big errors downstream



Part 3: Preparing files for analysis

Contamination


Contamination in your samples can lead to big errors downstream



Part 3: Preparing files for analysis

Contamination

Contamination in your samples can lead to big errors downstream



Outline

1. Different methods to acquire sequence data
2. Understand sequence file formats
3. Preparing files for analysis
 - **Tutorial looking at sequence data files and quality**

Tutorial:

Work through the tutorial associated with this session